

TRANSCEND™ with **PowerTraxx®**



The **Transcend™ Stair Chair with PowerTraxx®** provides maximum safety and power when transporting a patient up or down stairs.

Its electronically-controlled motor, direct drive transmission, and powered tracks push and pull patient loads of up to 550 pounds. The chair and track system carry the weight, increasing patient safety and reducing operator fatigue and risk of injury.

POWERCLIMB

with POWERTraxx technology can transport patients up and down an average of 40 flights of stairs with the touch of a button.

DESCENTCONTROL

limits track speed during descent or ascent to a safe, manageable level.

EZGLIDE®

provides a comfortable ride on almost any surface with 60% less shock than competitive stair chairs.

STAIRBRIGHT

offers integrated LED scene lights to illuminate surroundings up to 1 flight of stairs with two brightness options.



Motorized tracks

power the chair up and down stairs.



Extending footrest

provides added support and comfort for patients.



Rear lift handles

with improved ergonomics for lifting or carrying the chair.

**Preventative
Maintenance
Program**



*Preventative maintenance program - Parts, labor, travel, two annual preventative maintenance inspections, unscheduled service, and product health report. Replacement parts do not include: batteries,** accessories, other disposable or expendable parts, or damage due to misuse.*

**Batteries are included if accompanying Battery Plan is purchased.