

FERNO AIO Vacuum Mattresses & AS190 Splint Range

A range of vacuum products for effective immobilisation



FOR MORE DETAILED
INFORMATION SEE OVER

FERNO – is a world leading innovator and manufacturer of patient transport solutions.

Our product range is born of a closer understanding of our customers' needs and the practical demands of their work. Ferno's industry-leading **Vacuum Mattresses and Splints** typifies our passion for top quality, highly functional design... patient transfer products that don't just do the job, they do it **BETTER**.

Specifications

Product	Length (min)	Width	Height	Weight	Load Limit
AIO Mkl Mattress	2050mm	1150mm	55mm	8.7kg	220kg/34.6st/485lbs
AIO MkII Mattress	2050mm	1150mm	55mm	10.7kg	220kg/34.6st/485lbs
AS100 Arm Splint	550mm	350mm	30mm	1.55kg	N/A
AS120 Long Arm Splint	760mm	480/630mm	30mm	1.15kg	N/A
AS140 Leg Splint	11200mm	360/700mm	30mm	1.1kg	N/A
AS150 Long Leg Splint	1360mm	730mm	30mm	1.5kg	N/A
AS190 (Splint Kit)	800mm	210mm	310mm	4 kg	N/A

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk



FERNO AIO Vacuum Mattresses & AS190 Splint Range

Designed to provide a high degree of strength and stability, the Ferno range of splints & mattresses offer effective immobilisation.

- Highly flexible – moulds around the patient
- Flame retardant and easy to clean – fully welded without open seams or textiles
- Special valve design – no adapter needed
- Granulates disperse evenly to give the patient stable and optimal immobilisation



Part Numbers

038243600	All in One MkI 115cm Vacuum Mattress kit, c/w pump, repair kit & bag
27140400	AS192 Vacuum Mattress Kit c/w pump
271301001	All in One MkII 115cm Vacuum Mattress kit, c/w pump, repair kit & bag
081945900	AS190 Vacuum Splint Kit (AS100, AS120, AS140 Pump & Repair Kit)
090544500	AS100 Short Arm Vacuum Splint (incl. repair kit)
090544600	AS120 Long Arm Vacuum Splint (incl. repair kit)
090544700	AS140 Leg Vacuum Splint (incl. repair kit)
038266100	AS150 Extra Long Leg Vacuum Splint (incl. repair kit)

Accessories & Replacement parts

031347000	Replacement Orange bag for AS Vacuum Splint set
031331400	Replacement Carry Bag for Orange AIO and AS192 Vacuum Mattresses
083196600	Vacuum Splint repair Kit for Orange AIO Vacuum Mattress
271000701	Vacuum Pump with footrest for AIO & EasyFix Vacuum Mattress



FERNO - Vacuum Mattress & AS190 Splint Range. FPC003 10.15

call +44 (0) 1274 851 999
 email sales@ferno.co.uk
 www ferno.co.uk



Ferno (UK) Limited, Stubbs Beck Lane, Cleckheaton, West Yorkshire, BD19 4TZ England.