

FERNO EasyFix PLUS

Vacuum Mattress

Upgraded design for improved product durability and providing further enhanced immobilisation, removing the need to use a cervical collar



FOR MORE DETAILED
INFORMATION SEE OVER

Specifications

Length	Width (max)	Thickness	Weight	Load Limit
2080 mm	720 mm 1300 mm	45mm	8.2 kg (incl straps)	250 kg/39st/550lbs

FERNO – is a world leading innovator and manufacturer of patient transportation solutions.

Our product range is born of a closer understanding of our customers' needs and the practical demands of their work. Ferno's industry-leading **EasyFix PLUS Vacuum Mattress** typifies our passion for top quality, highly functional design... patient transfer products that don't just do the job, they do it **BETTER**.

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk

 **FERNO**[®]

FERNO EasyFix PLUS Vacuum Mattress

Following extensive user feedback we have upgraded and enhanced the EasyFix* to deliver even better head, neck and upper back immobilisation and enhanced durability of product when using in harsh operating situations.

FERNO's EasyFix PLUS provides an effective and secure vacuumed environment for a safer transfer of patients and casualties.

Designed to give **better product protection in unforgiving operating environments**, the ground-facing base surface is made from a durable low-friction material, with easier manoeuvrability and better patient handling.

A new V-shaped inner chamber system gives rigid and effective **C-spine immobilisation without the need to use a cervical collar**, which can cause extra trauma to the patient during positioning of a collar.

Designed for use in both rescue and medical emergency situations, the EasyFix PLUS vacuum mattress raises the level in patient immobilisation & transportation equipment, allowing a patient to be moved from scene all the way through to scanner without the need to transfer between different devices, in all environments and situations.

- **Improved design** – in the head and neck area for superior immobilisation of patient, with no collar required
- **Durable and highly effective** – base fibreglass material reduces friction making surface movement of the vacuum mattress easier
- **Integrated immobilisation** – integral head straps designed to remain in place (see images right)
- **High carrying capacity** – load limit 250kg
- **X-Ray and MRI applicable** – no need to move the patient to another device for diagnostics
- **Easy storage** – folds easily to take up 20% less space than other designs of vacuum mattress



*Standard FERNO EasyFix also available

Part Nos:

271511001	EasyFix PLUS Vacuum Mattress "Ready2Go kit" (incl. pump, repair kit & carry bag)
271512001	EasyFix PLUS Vacuum Mattress "BackPack Ready2Go kit" (incl. pump, repair kit & backpack)
271511000	EasyFix PLUS Vacuum Mattress (incl. restraints only)
271000708	Backpack for EasyFix PLUS, Red

FERNO - EasyFix PLUS Vacuum Mattress FPC002 11.17

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk



Ferno (UK) Limited, Stubs Beck Lane, Cleckheaton, West Yorkshire, BD19 4TZ England.