FERNO Models 11-B & **11-T** First Call Stretchers

Make removal down stairs and bariatric retrievals much simpler...with these innovative variations of our Model 11 first call stretcher

FERNO – is a world leading innovator and manufacturer of patient transportation and funeral equipment.
Our product range is born of a closer understanding of our customers' needs and the practical demands of their work.
Ferno's industry-leading Model 11 Stretcher typifies our passion for top quality, highly functional design...
funeral and mortuary products that don't just do the job, they do it BETTER.

Specifications

FOR MORE DETAILED

INFORMATION SEE OVER

Model 11-T Length Open	Length Folded	Width	Height Open	Depth Folded	Weight	Load Limit
1915 mm	920 mm	490 mm	180 mm	215 mm	10.5 kg	159 kg/25st/350lbs
Model 11-B Length Open	Length Folded	Width	Height Open	Depth Folded	Weight	Load Limit
1915 mm	975 mm	600 mm	210 mm	95 mm	7 kg	250 kg/39st/546lbs

call +44 (0) 1274 851 999 email sales@ferno.co.uk www ferno.co.uk



FERNO Models 11-B & **11-T** First Call Stretchers

Stay in complete control... with Bariatric and Tracked versions of our industry-leading First Call Stretcher

Model 11-B

A Bariatric First Call Stretcher for safe and dignified handling

- Exceptional ease of use lightweight but sturdy bariatric stretcher, with a load limit of 250kg (39 stone)
- Convenient for carrying or storage folds in half for storage
- Fast unfolding supplied with quick-release retaining straps
- Durable and hygienic polymer surface impervious to fluids
- Footrest included supplied as standard

Model 11-T

Perform removals down stairs with far greater confidence and control

Combines the innovative features of our Model 11 stretchers with a track system for ease of manoeuvring downstairs, minimising manual handling injury risk. SWL up to 159 kg (25 stone).

Part Nos:

60-0104-016/BLK	Model 11-T Stretcher with Track, incl. footplate, Black, 159 kg SWL
60-0104-010/BLK	Model 11-B Bariatric l Stretcher, incl. footplate, Black, 250 kg SWL

Options and Accessories:

031356800	Model 11-B Carry Bag, Black
031356800	Model 11-B Vinyl cover, Black
031346300	Model 11-T Carry Bag, Black
031356100	Model 11-T Box Style Carry Bag, Black

call +44 (0) 1274 851 999 **email** sales@ferno.co.uk **www** ferno.co.uk

Ferno (UK) Limited, Stubs Beck Lane, Cleckheaton, West Yorkshire, BD19 4TZ England.

FERNO Model 11 B & T FPC005 01.21

