

FERNO Paraguard 4-Point Lifting Sling

Designed and tested for lifting a casualty using Paraguard stretcher



The Paraguard 4-point lifting sling is designed for use with Paraguard rescue stretcher, tested and certificated for horizontal lifts up to the maximum SWL of Paraguard stretcher, load tested and rated up to 226 kgs. Incorporating metal carabiners that quickly attach to the 4 lifting points on the Paraguard stretcher, the Paraguard 4-point lifting sling is easy to deploy, giving Paraguard additional capability in various environments.

010773001	Paraguard Four Point Sling, for lifting
Weight	0.7 kg
Max SWL	226 kg
Lift-strap length	Head-end 80 cm, foot-end 100 cm

FERNO - Paraguard 4-point lifting sling FPC002 06.21

call +44 (0) 1274 851 999
email sales@ferno.co.uk
www ferno.co.uk

